

ANIKIT

Product Introduction

Simple and easy 1 minute DIY
check-up kit with smartphone



01. OVERVIEW


Empower Your Pet's Wellness Journey: Unlocking the potential to Prevent and Manage Health in Advance.



With the intention of providing even better care

Our Non-verbal companion animals

Due to the nuanced nature of non-verbal cues from our companions, it can be challenging to detect and prevent potential risks early and administer timely treatment



Understanding that frequent health check-ups might not always be possible

Since it's not easy going to hospital every time

Pet owners face limitations in conducting regular health check-ups for their companions due to cost and time constraints, also pets that fear clinics make the process even more challenging



Knowing that medical expenses can be burdensome

Increase in Consultation Costs

According to a press release by the Korea Consumer Federation, as of 2021, the average cost for a single medical consultation in 8,400 KRW which marks a 9.6% increase compared to 2019




Since pet owners are directly involved in managing their health

Take it into your own hands

Since silent-progressing illnesses don't cause immediate discomfort to pets, it's crucial for pet owners to take regular, proactive management measures

Empower Your Pet's Wellness Journey: Unlocking the potential to Prevent and Manage Health in Advance.

EXPENDITURE ON PET TREATMENT COST OVER PAST 2 YEARS



Therapeutic Expenses for Companion Animals



According to the '2021 South Korea Pet Report' released by KB Financial Group's Institute of Management Research, 71% of pet-owning households have incurred expenses for pet treatment costs over the past 2 years. Among these, dog owners spent an average of \$50 per visit, while cat owners paid an average of \$65 per visit



Silently Approaching Urinary Tract Disorders

Kidney deterioration goes unnoticed until it's significantly advanced, often leading to late-stage detection. Notably urinary tract disorders rank as the second leading cause of death in cats. Chronic kidney disease in cats is particularly fatal due to its propensity to trigger a range of complicating factors

Source: Excerpt from 'Daily Pet Clinic,' a specialized journal for veterinarians. [2021 Pet Ownership and Care Status]

Source: Excerpt from Daily Vet [Kidney Disease Ranks as second Leading Cause of Cat Mortality]

01. OVERVIEW

ANIKIT is available to be of assistance during moments like these

It's Important to Conduct a Health Check for Your Beloved Companion

When you haven't been consistent with regular health check-ups or when your pet displays behaviors that differ from their usual self, it's time to prioritize their health check



Existing Medical Conditions or Post Surgeries

If your companion has already received a diagnosis or undergone surgery at the hospital, it's important to regularly monitor their health at home



The Unusual Color of the Urine

If you notice any unusual color in your companion's urine, it could indicate various factors such as liver issues, urinary tract infections, and more



A Journey of Companionship Between a Healthy Companion and a Caring Owner


Ensuring your Companion's Health with Ease and Affordability, Right at Home

Since they are unable to verbally express pain, our companion's health demands vigilant attention and meticulous observation from its guardian. Diseases in companion animals can manifest without major warning signs. Especially in cases where ailments are not visibly apparent, it's very challenging to detect the problem. This is why routine check-ups are crucial. ANIKIT eases the burden by enabling seamless, at-home, and hassle-free regular pet health checks



02. PRODUCT INTRODUCTION


Developed a Highly Reliable Product Validated by National Authorities



01

Validation of Accuracy
and Precision

Tested and validated by KCL (Korea Construction Life Environment Testing Institute) for device and environmental conditions, achieving an accuracy of 98% and precision of 96%



02

Certification of Medical
Device for Animals

ANIKIT has been granted approval as a medical device for animals by the Animals and Plant Quarantine Agency (Approval No.342-001)



03

Urinalysis Test Strips
Patent Application

ANIKIT's urine test strips have been patented under the name of 'Smart Urinalysis Strips for Pets

02. PRODUCT INTRODUCTION

Detect multiple diseases with just one test results

Diabetes and Kidney Conditions



pH, Protein, Glucose, Ketones, Urine Specific Gravity

In the case of diabetes, it checks for detected glucose and ketones, as well as measures urine specific gravity, uric acid, and protein levels, which can indicate kidney function impairment

Liver Conditions



Bilirubin, Urobilinogen

In the case of impaired liver function, it interprets urobilinogen and bilirubin, which can be collected in urine, together to confirm liver health



Additional verifying factors

Ascorbic acid

Factors that could be Influenced by vitamins are collectively interpreted. So even if you detect warning signals in the test results, there's no need to worry

Urinary Tract Conditions

Hematuria, Leukocytes, Nitrite

You can collectively assess the health of the urinary tract system, particularly crucial for cats, by checking levels of hematuria, leukocytes, and nitrite

Manageable Conditions

Liver and Gallbladder Disorders, renal failure, diabetes, Ketosis, Hemolytic Anemia, Urinary Stones, Urinary Tract Infections, Bladder Stones Cystitis, etc

02. PRODUCT INTRODUCTION

Promptly and effortlessly review and record test results




Detecting 11 Test factors

You can review the levels and detailed information of Hematuria, Nitrite, Bilirubin, Glucose, Urobilinogen, pH, Ketones, Urine Specific Gravity, Protein, Leukocytes, and Ascorbic Acid



Effortless Testing with Automatic Scan in 1 Minute

After collecting urine, apply it to the test strip. One minute later, the camera app within the AhaPet app is automatically activated, facilitating easy scanning of the test strip



Real-Time Examination Result Review and Logging

Swiftly review the test results within the app after scanning. Gain detailed insights into each testing factor and keep records by session to track and monitor your companion's health status

03. INSTRUCTIONS

Using ANIKIT is Easy for Everyone

AhaPet Application Download


Android



iOS




01



Register your pet on the App

After downloading the 'AhaPet' app, sign up for an account. Once registered, add your pet and select the pet you want to test in the 'test' tab for examination

02



Collecting Urine Sample

Collect your companion's urine. Open the urine test strip and dip it in the urine. Carefully review the testing instructions in the app


03



Wipe Excess Urine

For accurate test results, if there's excessive urine on the test strip, gently press the side of the test strip against the absorbent pad to reabsorb the excess urine

04



Scan the Test stripe After 1 Minute

Place the urine test strip accurately in the center of the color chart and follow the app's instructions to automatically scan and instantly view the test results

Affordable and Convenient At-Home Health Care Possible with ANIKIT

